Ohřívač parafínu Weelko Melty

Model: WK-E005

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Model: WK-E005

GENERAL WARNING

Ensure that the voltage is equivalent to the mains voltage.

Do not use this device near a sink or other container filled with water.

Do not immerse the device in water or other liquid.

Always check the device and cable before use if you is not noticeable signs of damage. Do not use the device if it is damaged.

Do not leave the device unattended when it is switched on.

This device is not intended for use in children, as well as people with psychological, sensory or mental disabilities.

Switch off the device immediately after use.

Do not touch the heating element.

Do not place any flammable objects inside the tray or near the heater.

Make sure the paraffin level is not below "MIN" and not higher than the "MAX".

EC DECLARATION OF COMPLIANCE

This product meets the following requirements:

1. Conditions set out in the Low Voltage Directive: 2006/95/EC
2. Conditions set out in the Electromagnetic Compatibility Directive: 2004/108/EC

RECYCLING

Whenever any of our electrical products or their presentation packaging bear this symbol, it means that these products should not be treated as conventional household waste in Europe. To ensure proper treatment of this waste, dispose of it in accordance with local laws or as required for the disposal of electrical equipment. This initiative will help preserve save natural resources and improve standards for environmental protection in regard to the treatment of electrical waste.

WELCOME

Thank you very much for purchase. Please read the instructions carefully before operating the instrument. You can also ask a professional fur guidance to make sure you are using the instrument correctly. We will accept no responsibility for any accidents that occur due to incorrect operation.

INSTRUCTIONS FOR USE

1. Unpack the paraffin heater and place it on a stable surface.
2. Remove the lid, expand the paraffin and place it in the heater along with the rubber grille that is used to lay your hands or feet on the bottom. We recommend removing the grid after the paraffin bath and using it again in the melted paraffin after use.
3. Plug the device into the socket, close the lid and set the temperature button to "HI" - melting paraffin.
4. As soon as paraffin melts, switch the button to "LO". After slight cooling, check the temperature of wax by hand to see if it is possible to start with spa.
5. Remove from the hands and feet, everything you wear (eg. Rings, etc.). Wash your hands and feet thoroughly from dirt.
6. If the paraffin temperature is suitable for the bath, immerse your hands in paraffin 4-6 times for 2-3 seconds. A strong paraffin coating remains on the hands.
7. Apply foil gloves to maintain heat and for better absorption of nutrients from the paraffin bath. Put on terry gloves for 15-20 minutes.
8. After several treatments, the amount of paraffin drops to "MIN". In this case please add more wax to the "MAX" level. Please make sure that you always keep the level between the "MIN" and "MAX" marks.
9. Do not use the applied paraffin again.

CLEANING

Removal of paraffin from the container before cleaning:

- turn on the "HI" button - melt the paraffin and let the paraffin melt slightly after the bottom, then you can remove the entire paraffin block easily from the container,
- before cleaning clean the container and allow it to cool down.

When cleaning the appliance, always make sure that it is completely fine and unplugged.

After 100 treatments, you will need to clean the heater and replace the paraffin refill completely: Remove the paraffin from the heater and carefully clean the appliance from the inside and outside with a soft, clean cloth.

SPECIFICATION

Melting point: 65 °C

Operating temperature: 55 °C

Capacity: 2700 ml

Power: 150 W

Voltage: 220-240 V

Frequency: 50-60 Hz